


Medical Grade Oil Free Compressor



Medical Grade Oil Free Compressor



- ▶ Current: 4.8A
- ▶ Power: 1200W
- ▶ Volume flow: 150L/min @ 0bar
- ▶ Working pressure range: 5 up to 9.5bar
- ▶ Gas tank capacity: 45L
- ▶ Noise: 65-68dB
- ▶ High efficiency air dryer system

Medical Grade Oil Free Compressor



- ▶ Current: 9.6A
- ▶ Power: 2400W
- ▶ Volume flow: 300L/min @ 0bar
- ▶ Working pressure range: 5 up to 9.5bar
- ▶ Gas tank capacity: 55L
- ▶ Noise: 65-68dB
- ▶ High efficiency air dryer system

Centralized Air Compressor System



- ▶ Current: 18A
- ▶ Power: 4400W
- ▶ Volume flow: 600L/min @ 0bar
- ▶ Working pressure range: 5 up to 9.5bar
- ▶ Gas tank capacity: 150 - 300L

Medical Grade Oil Free Compressor



- ▶ Electric current: 9.6A
- ▶ Frequency: 2520W
- ▶ Power: 315L/min @ 0bar
- ▶ Working pressure range: 5 up to 9.5bar
- ▶ Gas tank capacity: 90L
- ▶ Noise: 68-75dB
- ▶ Smart Control System
- ▶ Optional: High efficiency air dryer system



Office # 301, 15 Paris Street, Swefieh
Amman, Jordan



www.periodentec.com



<https://web.facebook.com/periodentec>



@periodentec



+962 790663131

Periodent
Technologies

